

Tecotherm Neo

Thermoregulation and Monitoring



Total Body Cooling for infants suffering from Hypoxic Ischaemic Encephalopathy (HIE)

Tecotherm Neo

Thermoregulation and Monitoring

Total Body Cooling is a therapy for infants suffering from Hypoxic Ischaemic Encephalopathy (HIE) after Birth Asphyxia. The original Tecotherm was used in the biggest Randomised Control Trial for cooling infants with HIE, which showed improved neurologic outcomes for the cooled group¹. Now building on from that success, the Tecotherm Neo takes Total Body Cooling to the next level.

The innovative servo controlled design has instant feedback monitoring the infant's temperature every 2 seconds and making minute changes to the fluid to ensure that the infant's temperature remains stable. Simple alarms are on hand should the temperature deviate more than 0.5°C from the set temperature. The large colour display shows the temperature clearly and colour changes alert nursing staff to any problem.

The Tecotherm Neo can be programmed for a complete treatment cycle, the user can set the target temperatures, duration and even the rate of re-warming / cooling making it easy and simple to set up and use. Changes can be made at any time and all changes, set and measured parameters are recorded on to a memory card every minute for later analysis.



Stand by screen is easy to see across the room.



Clinician can set the duration of treatment, the rate of cooling / warming as well as the target temperature.



A lightweight and portable device, the Tecotherm Neo is the gold standard in Total Body Cooling.

Modes of Operation

The Tecotherm Neo has 3 modes of operation:

Servo-Control Complete Treatment Mode

Clinician can set the duration of treatment, the rate of cooling / warming as well as the target temperature. The Tecotherm Neo will complete the whole cycle without interruption and then maintain the final temperature.

Servo-Control Mode (Constant Rectal Temperature)

Clinician can set the target temperature, the time to get to target temperature and the length of time to maintain the target temperature. The Tecotherm Neo will complete the whole cycle without interruption.

Constant Mattress Temperature Mode

No Servo Control, useful for research work and where control of the mattress needs to be stable, or where experienced staff are available to monitor the patient.

The Tecotherm Neo takes Total Body Cooling to the next level

Flexible

- Clinician can set and adjust parameters at the start and throughout the treatment
- Customise and store treatment profiles (up to nine can be stored)
- Compatible with reusable and disposable accessories
- Optional use of skin probes

Easy to Use


- Large colour screen is easy to see across the room
- Small lightweight portable device
- The device restarts from the same settings in the event of power failure
- Maintains normothermia after rewarming

Compatible


- Data storage for allowing simple analysis of temperature profiles
- Export individual patient records to USB



Example of patient data collected from Tecotherm Neo



Technical Description and Specifications

Options	Cooling & Warming	Ambient conditions	
Patient weight max.	<50kg	Operation / Treatment Ambient Temperatures	+ 5°C to + 27°C
Dimensions	375mm (width) x 190mm / 215mm (height) x 310mm (depth)	Operation / Treatment Relative Humidity	10% to 75%, not condensating
Weight without accessories	approx. 7.2kg	System safety	
Central cooling module	Thermoelectrically based module	Protection class	Class 1, Risk Class II b, Type BF 
Treatment temperature control ranges	Mattress + 12°C to + 39°C Rectal BCT + 30°C to + 38°C	Standards	DIN EN 60601-1 DIN EN 60601-1-2 DI 60601-1-6 DIN EN 60601-1-8 DIN EN 60601-1-10 DIN EN 60601-2-35 E/F
Hydraulic circulation system			
System pressure max.	0,5 bar		
Flow rate without / with mattress	500 ml / min (shorted) / up to 300 in use		
Internal fluid reservoir capacity	approx. 250 ml		
Circulating fluid	Sterile Water		
Connectors / Couplings	Quick Disconnect Couplings		
Fill up / Refill	Fill up set		
Electrical parameters			
Supply voltage / Mains	100-130V and 200-240V, 50-60 Hz		
Power consumption	max. 350 W		
Patient safety / alarms			
Lower temperature alarm limit	+10°C		
Upper temperature alarm limit	+41°C		
Set temperatures, lower limit	+12°C		
Set temperatures, upper limit	+39°C		

References

1. Azzopardi et al (2009) Moderate Hypothermia to Treat Perinatal Asphyxial Encephalopathy, N Engl J Med 2009; 361:1349-1358

The details given in this brochure are correct at time of going to press. The company reserves the right to improve the products shown.

500-126(1) 2018-11



Distributed by

Inspiration Healthcare Limited

2 Satellite Business Village, Fleming Way,
Crawley, West Sussex, RH10 9NE
United Kingdom

T +44 (0)1455 840 555

F +44 (0)1455 841 464

E info@inspiration-healthcare.com

W inspiration-healthcare.com

Autorizované zastoupení pro ČR:

DARTiN

DARTIN, spol. s r.o.
Suchdolská 688, 252 62 Horoměřice
Tel./fax: 241 470 361, 261 216 076
www.dartin.cz
e-mail: info@dartin.cz

Autorizované zastúpenie pre SR:

DARTiN

DARTIN Slovensko, s.r.o.
Ľ. Podjavorinskej 25, 080 05 Prešov
Tel./fax: 051/7734403, 7723634
www.dartin.sk
e-mail: info@dartin.sk