

DATASHEET

Eye protector biliSAVE[®]


Page | 1

| | |
|------------------------|--|
| Art.-No.: | NA421644/1 • NA421643/1 • NA421642/1 |
| Art.-Description: | Eye protector biliSAVE [®] for phototherapy |
| Manufacturer: | Kanmed AB |
| Exclusive Distributer: | MedCare Visions [®] GmbH • Franz-Lehner-Str. 3 • 85716 Unterschleißheim |

Delivery scope

| | | |
|------------|---|---------|
| NA421644/1 | biliSAVE [®] size S - for head circumference 20 - 28cm | 50 pcs. |
| NA421643/1 | biliSAVE [®] size M - for head circumference 24 - 33cm | 50 pcs. |
| NA421642/1 | biliSAVE [®] size L - for head circumference 30 - 38cm | 50 pcs. |

Product picture and quick guide



Technical data and features

| | |
|-----------------------|---|
| Material texture: | Spun laced non-woven |
| Product components: | Nylon (Velcro) |
| Application type: | Single patient use |
| Product construction: | Eye protector with Velcro closure |
| Intended Use: | The biliSAVE [®] neonatal eye protector is intended to be used during phototherapy of new-born babies to prevent retinal damage and to provide optimal irradiated surface. It shall be used by adult users with normal cognition and in a professional care or home care environments. |
| Warning: | Single patients use only. Reuse for other patients may cause cross infection. The device is intended for 24h of use, prolonged use may cause skin irritation. |

Autorizované zastoupení pro ČR:

DARTIN

DARTIN, spol. s r.o.
Suchdolská 688, 252 62 Horoměřice
Tel./fax: 241 470 361, 261 216 076
www.dartin.cz
e-mail: info@dartin.cz