

Golvo™ 9000/LowBase™ mobile lift

Designed for versatility





Golvo 9000 LowBase

The Golvo 9000 lift offers the functional advantages of an overhead lift with the flexibility of a mobile lift. In combination with the extensive range of Liko accessories, Golvo is ideal for common lifting situations across a variety of care environments.

Features and benefits

Vertical lifting motion with flexible lift strap is ideal for horizontal lifting, ambulation and lifting from the floor

Parallel base-width adjustment contributes to easier lifting of patients or residents to/from chairs, beds and the floor

LowBase version allows for even greater clearance under low beds or stretchers versus the standard model

Armrests help to **support patients during ambulation**, as well as making it easier for caregivers to manoeuvre the lift

User display highlights device usage, including lift count, service required indicator and overloads

Optional **quick-release hook system** means caregivers can easily and rapidly change a sling bar or install a scale

Ergonomic hand control for lift operation displays both battery usage and overload indicators

An extensive assortment of accessories is available including slings, sling bars and scales


Two lifting speeds accommodate a variety of different lifting situations

Golvo 9000/LowBase mobile lift









Lift arm LowBase - Base height

Hand control

Technical Specifications	
Maximum Load	200 kg
Weight - 9000 with lead battery with Li lon battery	41.6 kg 40.2 kg
Weight - 9000 LowBase with lead battery with Li lon battery	45 kg 43.6 kg
Height 9000 LowBase	2090 mm 2090 mm
Protection class	IPX4
Hand control protection class	IPX4
Slingbar options standard optional	Universal 450 slingbar (450mm) Mini slingbar 220, Universal slingbar 350, Universal TwinBar 670, Crossbars
Base - 9000 height length width	105 mm 1185 mm 735 - 1020 mm

Base - 9000 LowBase	
height	60 / 105 mm
length	1185 mm
width	745 - 1028 mm
Base - width expansion	electric
Castors	
front - 9000	75 mm double band
front - LowBase	46 mm double band
rear	75 mm double band
Material	aluminium
Emergency lifting	electric
Emergency lowering	manual and electric
Lifting speed (without load)	27 mm / sec or
	37 mm / sec
Lifting interval	1245 mm
Maximum lifting height	1816 mm
Sound level	52.8 dB(A)
Battery option	Lead acid Gel or
	Lithium Ion
Charger	inbuilt
Hand controls	wired
Product number	
9000	2000045
LowBase	2000049

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health. Around the world the people, programs, and products of Hill-Rom work towards one mission: Enhancing outcomes for patients and their caregivers.

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

Golvo 9000 mobile lift is intended to be used for the lifting and transfer of patients in a variety of care settings.

Class I

Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Autorizované zastoupení pro ČR:



DARTIN, spol. s r.o. Suchdolská 688, 252 62 Horoměřice Tel./fax: 241 470 361, 261 216 076 www.dartin.cz e-mail: info@dartin.cz Autorizované zastúpenie pre SR:



DARTIN Slovensko, s.r.o. L'. Podjavorinskej 25, 080 05 Prešov Tel./fax: 051/7734403, 7723634 www.dartin.sk e-mail: info@dartin.sk

Enhancing outcomes for patients and their caregivers:

